

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB17740 | Gene ID: | 57532 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 76kDa | Conjugate: | Unconjugated |
| Calculated MW: | 76kDa | Isotype: | IgG |

Immunogen Information

Background: Enables RNA binding activity. Located in cytoplasmic stress granule; cytosol; and nuclear body. Part of polysomal ribosome.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | PIG1; NUFP2; 82-FIP; FIP-82; 182-FIP; NUFIP2 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 456-695 of human NUFIP2 (NP_065823.1). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |