

NCEH1 Rabbit Polyclonal Antibody



CAB17741

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

46kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57552

Uniprot

Q6PIU2

Synonyms:

AADACL1; NCEH; NCEH1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 90-440 of human NCEH1 (NP_065843.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - NCEH1 Polyclonal Antibody (CAB17741)