

Product Information

Product SKU:	CAB17747	Gene ID:	60493	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	87kDa	Conjugate:	Unconjugated
Calculated MW:	87kDa	Isotype:	IgG

Immunogen Information

Background: Enables rRNA binding activity. Involved in mitochondrial RNA processing. Located in mitochondrial nucleoid and ribonucleoprotein granule.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	dJ1187M17.5; FASTKD5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 495-764 of human FASTKD5 (NP_068598.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.