

# FASTKD5 Rabbit Polyclonal Antibody



CAB17747

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

87kDa

### Calculated MW:

### Applications:

WB IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

60493

#### Uniprot

Q7L8L6

#### Synonyms:

dJ1187M17.5; FASTKD5

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 495-764 of human FASTKD5 (NP\_068598.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - FASTKD5 Rabbit pAb (CAB17747)