

Product Information

Product SKU:	CAB17751	Gene ID:	63874	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	39kDa	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable lysophosphatidic acid acyltransferase activity and lysophospholipase activity. Predicted to be involved in N-acylphosphatidylethanolamine metabolic process; lipid homeostasis; and phosphatidic acid biosynthetic process. Predicted to act upstream of or within N-acylethanolamine metabolic process. Predicted to be located in endoplasmic reticulum membrane. Predicted to be active in lipid droplet and mitochondrion.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: ABH4; ABHD4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 223-342 of human ABHD4 (NP_071343.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.