

AZI2 Rabbit Polyclonal Antibody



CAB17756

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

26kDa/28kDa/44kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

AZI2, or NAP1, contributes to the activation of NFKB (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK (TBK1; MIM 604834) (Fujita et al., 2003 [PubMed 14560022]).[supplied by OMIM, Mar 2008]

Immunogen information

Gene ID:

64343

Uniprot

Q9H6S1

Synonyms:

AZ2; NAP1; TILP; AZI2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic peptide of Human AZI2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
