

# FNDC4 Rabbit Polyclonal Antibody



CAB17758

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

25kDa

### Applications:

IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

64838

#### Uniprot

Q9H6D8

#### Synonyms:

FNDC4; FRCP1

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 45-167 of human FNDC4 (NP\_073734.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry - FND4 Rabbit pAb (CAB17758)