

# PRKRIP1 Rabbit Polyclonal Antibody



CAB17769

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

21kDa

### Calculated MW:

### Applications:

WB IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

79706

#### Uniprot

Q9H875

#### Synonyms:

C114; KRBOX3; PRKRIP1

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 67-134 of human PRKRIP1 (NP\_078929.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - PRKRIP1 Polyclonal Antibody (CAB17769)

