
Product Information

Product SKU:	CAB17770	Gene ID:	79720	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	31kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background: Enables calcium-dependent protein binding activity. Involved in positive regulation of viral budding via host ESCRT complex. Located in endosome membrane; midbody; and plasma membrane. Part of ESCRT I complex.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: VPS37B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 60-170 of human VPS37B (NP_078943.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.