

VPS37B Rabbit Polyclonal Antibody



CAB17770

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

31kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79720

Uniprot

Q9H9H4

Synonyms:

VPS37B


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 60-170 of human VPS37B (NP_078943.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - VPS37B Polyclonal Antibody (CAB17770)