

Product Information

|              |          |               |        |             |             |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB17774 | Gene ID:      | 79954  | Size:       | 20uL, 100uL |
| Clone No:    | -        | Host Species: | Rabbit | Reactivity: | Human,Mouse |

Additional Information

|                |       |            |              |
|----------------|-------|------------|--------------|
| Observed MW:   | 90kDa | Conjugate: | Unconjugated |
| Calculated MW: | 80kDa | Isotype:   | IgG          |

Immunogen Information

**Background:** Enables RNA binding activity. Predicted to be involved in maturation of SSU-rRNA from tricistronic rRNA transcript (SSU-rRNA, 5.8S rRNA, LSU-rRNA). Located in nucleolus.

|                       |   |
|-----------------------|---|
| Recommended Dilution: | WB,1:500 - 1:2000   |
| Synonyms:             | PQBP5; NOL10  |
| Purification Method:  | Affinity purification   |
| Immunogen:            | Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human NOL10 (NP_079170.2). |
| Storage:              | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.                    |