

Product Information

Product SKU:	CAB17776	Gene ID:	80067	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	59kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a nuclear transmembrane protein that associates with cullin 4A/damaged DNA binding protein 1 ubiquitin ligase complex. Mutations in this gene are associated with Woodhouse-Sakati syndrome. Alternate splicing results in multiple transcript variants.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	C2orf37; C20orf37; DCAF17
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 235-335 of human DCAF17 (NP_079276.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.