

# PUS7L Rabbit Polyclonal Antibody



CAB17783

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

81KDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

83448

#### Uniprot

Q9H0K6

#### Synonyms:

PUS7L

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PUS7L (NP\_112582.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - PUS7L Rabbit pAb (CAB17783)