

COG3 Rabbit Polyclonal Antibody



CAB17785

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

94kDa

Calculated MW:

50kDa/94kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a component of the conserved oligomeric Golgi (COG) complex which is composed of eight different subunits and is required for normal Golgi morphology and localization. Defects in the COG complex result in multiple deficiencies in protein glycosylation. The protein encoded by this gene is involved in ER-Golgi transport.[provided by RefSeq, Jun 2011]

Immunogen information

Gene ID:

83548

Uniprot

Q96JB2

Synonyms:

COG3; SEC34

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human COG3 (NP_113619.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - COG3 Polyclonal Antibody (CAB17785)