
Product Information

Product SKU:	CAB17785	Gene ID:	83548	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	94kDa	Conjugate:	Unconjugated
Calculated MW:	94kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a component of the conserved oligomeric Golgi (COG) complex which is composed of eight different subunits and is required for normal Golgi morphology and localization. Defects in the COG complex result in multiple deficiencies in protein glycosylation. The protein encoded by this gene is involved in ER-Golgi transport.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: SEC34; COG3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human COG3 (NP_113619.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.