

SLC25A18 Rabbit Polyclonal Antibody



CAB17786

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

33kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

83733

Uniprot

Q9H1K4

Synonyms:

GC2; SLC25A18

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 90-220 of human SLC25A18 (NP_113669.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
