

MARVELD1 Rabbit Polyclonal Antibody



CAB17787

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

19kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

83742

Uniprot

Q9BSK0

Synonyms:

GB14; MARVD1; MRVLDC1; bA548K23.8; MARVELD1


Immunogen:

A synthetic peptide of Human MARVELD1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MARVELD1 Polyclonal Antibody (CAB17787)