

Product Information

Product SKU:	CAB17787	Gene ID:	83742	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	19kDa	Conjugate:	Unconjugated
Calculated MW:	19kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be a structural constituent of myelin sheath. Predicted to be involved in myelination. Predicted to be located in several cellular components, including cytoskeleton; nucleus; and plasma membrane. Predicted to be integral component of membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	GB14; MARVD1; MRVLDC1; bA548K23.8; MARVELD1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-50 of human MARVELD1 (NP_113672.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.