

CAB17790

Product Information

Product SKU:	CAB17790	Gene ID:	84258	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	63kDa	Conjugate:	Unconjugated
Calculated MW:	63kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable several functions, including phospholipid binding activity; protein dimerization activity; and syntaxin binding activity. Involved in positive regulation of dendrite extension. Located in endosome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SytIII; SYT3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 449-590 of human SYT3 (NP_115674.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.