

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB17795 | Gene ID: | 84654 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Rat |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 50kDa | Conjugate: | Unconjugated |
| Calculated MW: | 49kDa | Isotype: | IgG |

Immunogen Information

Background: This gene encodes a bHLH-zip transcription factor which functions in the mitogen-activate protein kinase (MAPK) signaling pathway. Because of its role in the upregulation of cell proliferation and tumorigenesis, this gene may serve as a target for Ras-induced tumor treatments.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:1000 |
| Synonyms: | CILD38; NYD-TSP1; PPP1R148; SPZ1 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-430 of human SPZ1 (NP_115956.3). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |