
Product Information

Product SKU:	CAB17797	Gene ID:	84969	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: Enables transcription coactivator activity. Involved in positive regulation of transcription by RNA polymerase II. Located in nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: GCX1; GCX-1; C20orf100; dJ495O3.1; dJ1108D11.2; TOX2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human TOX2 (NP_116272.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.