

# TOX2 Rabbit Polyclonal Antibody



CAB17797

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

70kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

84969

#### Uniprot

Q96NM4

#### Synonyms:

C20orf100; GCX-1; GCX1; dJ1108D11.2; dJ495O3.1; TOX2

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human TOX2 (NP\_116272.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - TOX2 Polyclonal Antibody (CAB17797)