

Product Information

Product SKU:	CAB17802	Gene ID:	90634	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	29kDa	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

Background: N4BP2L1 (NEDD4 Binding Protein 2 Like 1) is a Protein Coding gene. Diseases associated with N4BP2L1 include Inflammatory Leiomyosarcoma. Gene Ontology (GO) annotations related to this gene include kinase activity. An important paralog of this gene is N4BP2.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CG018; N4BP2L1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human N4BP2L1 (NP_001073159).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.