

# TRIM47 Rabbit Polyclonal Antibody



CAB17803

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

70kDa

### Calculated MW:

### Applications:

WB IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

91107

#### Uniprot

Q96LD4

#### Synonyms:

GOA; RNF100; TRIM47

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 356-530 of human TRIM47 (NP\_258411.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - TRIM47 Polyclonal Antibody (CAB17803)