

Product Information

Product SKU:	CAB17803	Gene ID:	91107	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin-protein transferase activity. Involved in protein ubiquitination. Located in cytosol.

Recommended Dilution:	WB,1:100 - 1:500 IF/ICC,1:50 - 1:200
Synonyms:	GOA; RNF100; TRIM47
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 356-530 of human TRIM47 (NP_258411.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.