

DEPDC7 Rabbit Polyclonal Antibody



CAB17805

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

58kDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

91614

Uniprot

Q96QD5

Synonyms:

TR2; dJ85M6.4; DEPDC7


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human DEPDC7 (NP_001070710.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - DEPDC7 Polyclonal Antibody (CAB17805)