

Product Information

Product SKU:	CAB17805	Gene ID:	91614	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	58kDa	Conjugate:	Unconjugated
Calculated MW:	58kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in intracellular signal transduction.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	TR2; dJ85M6.4; DEPDC7
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human DEPDC7 (NP_001070710.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.