

# PAXBP1 Rabbit Polyclonal Antibody



CAB17812

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse

## Protein Background

This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009]

## Immunogen information

### Gene ID:

94104

### Uniprot

Q9Y5B6

### Synonyms:

BM020; C21orf66; FSAP105; GCFC; GCFC1; PAXBP1

### Immunogen:

A synthetic peptide of Human PAXBP1.

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---