

---

## Product Information

<b>Product SKU:</b>	CAB17815	<b>Gene ID:</b>	114791	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Rat

---

## Additional Information

<b>Observed MW:</b>	118kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	118kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables microtubule binding activity. Involved in microtubule nucleation. Located in centrosome and cytosol. Part of gamma-tubulin large complex.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** GCP5; TUBGCP5

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TUBGCP5 (NP\_443135.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.