

CAB17818

Product Information

Product SKU:	CAB17818	Gene ID:	123811	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	20kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables identical protein binding activity. Involved in cilium assembly. Located in centriolar satellite and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FOPNL; FOR20; C16orf63; PHSECRG2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-100 of human FOPNL (NP_653201.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.