
Product Information

Product SKU:	CAB17820	Gene ID:	125058	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	86kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa	Isotype:	IgG

Immunogen Information

Background: Enables GTPase activator activity. Involved in regulation of receptor recycling. Located in cytosol and early endosome.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: TBC1D16

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TBC1D16 (NP_061893.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.