

# TBC1D16 Rabbit Polyclonal Antibody



CAB17821

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

IHC IF

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:100 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

125058

#### Uniprot

Q8TBP0

#### Synonyms:

TBC1D16

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TBC1D16 (NP\_061893.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---