

Product Information

Product SKU:	CAB17823	Gene ID:	131177	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

Background: Involved in negative regulation of insulin secretion. Predicted to be located in extracellular region. Predicted to be integral component of membrane. Predicted to be active in extracellular space.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	EF7; OIT1; FAM3D
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 24-204 of human FAM3D (NP_620160.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.