

# FAM3D Rabbit Polyclonal Antibody



CAB17823

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

25kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

131177

#### Uniprot

Q96BQ1

#### Synonyms:

EF7; OIT1; FAM3D

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 24-204 of human FAM3D (NP\_620160.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - FAM3D Polyclonal Antibody (CAB17823)

