

CCDC12 Rabbit Polyclonal Antibody



CAB17835

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

19kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

151903

Uniprot

Q8WUD4

Synonyms:

CCDC12

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 44-123 of human CCDC12 (NP_653317.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
