

CAMSAP1 Rabbit Polyclonal Antibody



CAB17839

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

240kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

157922

Uniprot

Q5T5Y3

Synonyms:

CAMSAP1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 350-650 of human CAMSAP1 (NP_056262.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - CAMSAP1 Rabbit pAb (CAB17839)