

APLF Rabbit Polyclonal Antibody



CAB17842

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

57kDa

Calculated MW:

56kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

C2ORF13 is a component of the cellular response to chromosomal DNA single- and double-strand breaks (Iles et al., 2007 [PubMed 17353262]).[supplied by OMIM, Mar 2008]

Immunogen information

Gene ID:

200558

Uniprot

Q8IW19

Synonyms:

APFL; C2orf13; PALF; Xip1; ZCCHH1; APLF


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human APLF (NP_775816.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - APLF Rabbit pAb (CAB17842)