

Product Information

Product SKU:	CAB17843	Gene ID:	202052	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable Hsp70 protein binding activity. Predicted to be involved in cellular response to misfolded protein; chaperone cofactor-dependent protein refolding; and ubiquitin-dependent ERAD pathway. Predicted to be integral component of membrane. Predicted to be active in endoplasmic reticulum membrane.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	DNAJC18
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 260-350 of human DNAJC18 (NP_689899.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.