

Cyclin Y Rabbit Polyclonal Antibody



CAB17845

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

39kDa

Calculated MW:

33kDa/36kDa/39kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Cyclins, such as CCNY, control cell division cycles and regulate cyclin-dependent kinases (e.g., CDC2; MIM 116940) (Li et al., 2009 [PubMed 18060517]).[supplied by OMIM, May 2009]

Immunogen information

Gene ID:

219771

Uniprot

Q8ND76

Synonyms:

C10orf9; CBCP1; CCNX; CFP1; CCNY; cyclin-Y

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 60-120 of human CCNY (NP_859049.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - Cyclin Y Rabbit pAb (CAB17845)

