

Product Information

Product SKU:	CAB17846	Gene ID:	256364	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	95kDa	Conjugate:	Unconjugated
Calculated MW:	95kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable microtubule binding activity. Involved in mitotic metaphase plate congression and regulation of mitotic spindle assembly. Located in several cellular components, including midbody; mitotic spindle microtubule; and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ELP95; EMAP3; EMAP95; EML3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 300-530 of human EML3 (NP_694997.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.