

ARHGAP30 Rabbit Polyclonal Antibody



CAB17847

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

119kDa

Calculated MW:

Applications:

WB IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

257106

Uniprot

Q7Z6I6

Synonyms:

ARHGAP30


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 485-585 of human ARHGAP30 (NP_859071.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ARHGAP30 Rabbit pAb (CAB17847)