

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-------------|
| Product SKU: | CAB17847 | Gene ID: | 257106 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse |

Additional Information

| | | | |
|----------------|--------|------------|--------------|
| Observed MW: | 119kDa | Conjugate: | Unconjugated |
| Calculated MW: | 119kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable GTPase activator activity. Predicted to be involved in small GTPase mediated signal transduction. Located in intracellular membrane-bounded organelle.

| | |
|-----------------------|--|
| Recommended Dilution: | WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 |
| Synonyms: | ARHGAP30 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 485-585 of human ARHGAP30 (NP_859071.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |