
Product Information

Product SKU:	CAB17853	Gene ID:	320460	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	23kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: Acts upstream of or within negative regulation of BMP signaling pathway and positive regulation of neuron differentiation. Located in extracellular space. Part of AMPA glutamate receptor complex. Is expressed in brain; central nervous system; hippocampus; hypothalamus; and thalamus. Orthologous to human VWC2L (von Willebrand factor C domain containing 2 like).

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: Brl; A830006F12Rik; VWC2L

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-50 of mouse VWC2L (NP_796138.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.