
Product Information

Product SKU:	CAB17862	Gene ID:	219333	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Enables cysteine-type endopeptidase activity and thiol-dependent deubiquitinase. Involved in protein deubiquitination. Predicted to be located in nucleoplasm. Predicted to be active in cytosol and nucleus.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: UBH1; USP12L1; USP12

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP12 (NP_872294.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.