

sus1 Rabbit Polyclonal Antibody



CAB17870

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

93kDa

Applications:

WB

Reactivity:

Zea mays

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

542247

Uniprot

P49036

Synonyms:

ss2; GRMZM2G152908; sus1

Immunogen:

A synthetic peptide of Zea mays sus1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
