

# ZC3H12D Rabbit Polyclonal Antibody



CAB17881

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

70kDa

### Calculated MW:

25kDa/36kDa/58kDa

### Applications:

WB IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

340152

#### Uniprot

A2A288

#### Synonyms:

C6orf95; MCPIP4; TFL; dJ281H8.1; p34; ZC3H12D

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 150-230 of human ZC3H12D (NP\_997243.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - ZC3H12D Rabbit pAb (CAB17881)