

Product Information

Product SKU:	CAB17881	Gene ID:	340152	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	58kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable endoribonuclease activity and mRNA binding activity. Involved in negative regulation of G1/S transition of mitotic cell cycle and negative regulation of cell growth. Located in P-body and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	TFL; p34; MCPIP4; C6orf95; dJ281H8.1; ZC3H12D
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 150-230 of human ZC3H12D (NP_997243.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.