

Pan PKC Rabbit Polyclonal Antibody



CAB17921

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

83kDa

Calculated MW:

46kDa/56-68kDa/74-83kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5578/5579/5580/5581/5582/5583/5584/5588/5590

Uniprot

P17252/P05771/Q05655/Q02156/P05129/P24723/P41743/Q04759/Q05513

Synonyms:

Immunogen:

A synthetic peptide of human Pan PKC.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - Pan PKC Rabbit pAb (CAB17921)

