

PXDN Rabbit Polyclonal Antibody



CAB17929

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

165KDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a heme-containing peroxidase that is secreted into the extracellular matrix. It is involved in extracellular matrix formation, and may function in the physiological and pathological fibrogenic response in fibrotic kidney. Mutations in this gene cause corneal opacification and other ocular anomalies, and also microphthalmia and anterior segment dysgenesis.

Immunogen information

Gene ID:

7837

Uniprot

Q92626

Synonyms:

PXN; VPO; MG50; PRG2; ASGD7; COPOA; D2S448; D2S448E; PXDN


Immunogen:

Recombinant protein of human PXDN.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PXDN Rabbit pAb (CAB17929)