

# WRKY8 Rabbit Polyclonal Antibody



CAB17934

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

## Protein Background

### Immunogen information

#### Gene ID:

834678

#### Uniprot

### Applications:

WB

#### Synonyms:

WRKY8

### Reactivity:

Arabidopsis thaliana

#### Immunogen:

Recombinant protein of human WRKY8.

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Product Images

---